HR/BIH Korisnički priručnik

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EN

The device bears CE mark indicating its conformity with the provisions of the Council Directive 93/42/EEC concerning medical devices and fulfills the essential requirements of Annex I of this directive. The CT14 Continuous Glucose Monitoring System is manufactured by Zhejiang POCTech Co., Ltd. © Zhejiang POCTech Co., Ltd. All rights reserved. This User Manual may change without notice. When this happens, a new issue of the User Manual will be released.

Definitions of Safety Notes

 **Warning**

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in serious consequences.

 **Caution**

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the medical device or other property.

Note

A NOTE provides additional information intended to avoid inconvenience during operation.

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Contents

Section 1. Safety Information	88
1.1. Important Safety Information	88
1.2. Labels	88
Section 2. Product Introduction	89
2.1. Intended Use	89
2.2. Contraindications	89
2.3. System Introduction	89
2.4. System Components	89
2.5. Receiver Menu Description	91
Section 3. Preparation before Use	93
3.1. Preparation before Use	93
Section 4. Product Operation and Notice	94
4.1. Using a New Sensor	94
4.2. Initialize the Sensor	98
4.3. Sensor Initialization and First Ref BG	98
4.4. Events Input	99
4.5. Features of the Touch-Screen Display	100
4.6. Care and Maintenance during Monitoring	100
4.7. Ending Your Continuous Glucose Monitoring	101
4.8. Waste Disposal	101
Section 5. Troubleshooting and Maintenance	102
5.1. Troubleshooting(with Error Code)	102
5.2. General Maintenance of Device -Inspection and Maintenance	104
Section 6. Specifications and EMC Declaration	104
6.1. Specifications	104
6.2. EMC Declaration	105
6.3. EU Radio Equipment Directive (RED) Declaration	107
Section 7. Warranty	107
7.1. Contacts and Numbers	107
7.2. Warranty	107

Section 1. Safety Information

1.1. Important Safety Information

Read this entire manual before using your Continuous Glucose Monitoring System. If you do not understand something or have any questions, please ask your Diabetes Management Team or another health care provider. Contraindications, warnings, precautions, cautions, and other important safety information can be found in this User Manual.

Warning

This device is not designed to replace a blood glucose meter. The system must be used with a blood glucose meter in compliance with regulations.

- Treatment decisions should not be based solely on results from this device. You must confirm with a blood glucose meter before making therapeutic adjustments.
- Symptoms related to low or high blood glucose levels should not be ignored. If you have symptoms of low or high glucose, use your blood glucose meter to verify the results.
- You should update the device's calibration according to the instruction in this manual to ensure device performance. The performance of the device when calibrated less frequently than the recommendation has not been studied.
- The CGM product should be used under the supervision of a physician.

Caution

- Always wash hands with soap and water before opening the sensor package. After opening the package, avoid touching the adhesive area.
- Before inserting the sensor, always clean the skin at the sensor insertion location with a topical antimicrobial solution such as isopropyl alcohol. Do not apply the sensor until the cleaned area is dry.
- Establish a rotation schedule for choosing each new sensor location. Avoid sensor locations that are constrained by clothing, accessories, or subjected to rigorous movement during exercise.
- Avoid injecting insulin or placing an insulin pump infusion set within 3 inches of a sensor.
- The sensor is sterile in its unopened, undamaged package. Do not use any sensor if its sterile package has been previously damaged or opened.
- For calibration you must enter the exact finger stick reading that your blood glucose meter displays. Enter all finger stick readings for calibration within 5 minutes. Entering incorrect finger stick readings that occurred more than 5 minutes ago will affect device performance.

1.2. Labels:

Symbol label	Title of symbol
	Consult instructions for use
	Refer to instruction manual / booklet
STERILE R	Sterilized using irradiation
	Do not re-use
	Type BF applied part
	Direct current
	Non-ionizing electromagnetic radiation
	Symbol for the marking of electrical and electronic equipment
LOT	Batch code
	Use-by date
	Temperature limit (2 °C ~ 30 °C)
IP27	Degrees of protection against water indicated
	Manufacturer
	Do not use if package is damaged
	Date of manufacture
EC REP	Authorized representative in the European Community

Section 2. Product Introduction

2.1. Intended Use

The device is indicated for detecting trends and tracking patterns in adults with diabetes. The system is intended for use by patients. This device is not designed to replace a blood glucose meter. The device must be used with a blood glucose meter. Treatment decisions should not be based solely on results from the system. Treatment decisions must be confirmed with a blood glucose meter before making therapeutic adjustments.

Indications

The system is indicated for detecting trends and tracking patterns in adults with diabetes. It is intended for use by patients.

2.2. Contraindications

CT-Scan, X-Ray or MRI: If you must be examined with CT, X-ray, or MRI imaging, remove the transmitter and sensor. Start a new sensor after the examination.

2.3. System Introduction:

Principle of Operation

The sensor works through enzymatic electrochemical reactions with glucose. The enzyme, glucose oxidase, allows glucose to be converted into an electronic signal. The receiver converts the sensor signal to a glucose reading that is displayed to the user.

The CT14 allows you to continually see your readings every 3 minutes for up to 14 days.

Reagents

Glucose oxidase. The enzyme glucose oxidase is derived from a microorganism *Aspergillus niger* purified and dried according to Type VII-S guidelines. Direct tissue contact with the glucose oxidase is prevented by an outer sensor membrane layer, therefore the risk of tissue reactions from glucose oxidase is considered to be minimal.

2.4. System Components

The system consists of three parts (Fig.1), a glucose sensor, a transmitter, and a receiver.


Figure 1. System: Sensor, Transmitter, and Receiver

- The sensor is a four electrode flexible tiny strip assembled with an inserter, packed in a sterile package for direct use.


Figure 2. The Sensor Assembly (Sensor & Inserter)

Part		
1	1.1	CT-202 Sensor
2	2.1	CT-100C ₁₄ Transmitter
3	3.1	CT-100BD Receiver
4	4.1	User Manual
	4.2	Quick Start Guide
	4.3	Accessory: USB (Type C) Cable
	4.4	Accessory: Battery Lid
	4.5	Accessory: CR1620 Lithium ion button battery for CT-100C14 Transmitter(Optional)

• Additional Supplies Needed:

Users should have a blood glucose meter readily accessible to make therapeutic adjustments and calibrate the CGM.

Users can choose to receive data using the receiver or POCTech CGM App. Users can download the POCTech CGM App through the mobile application market. Used with the transmitter and sensor, the POCTech CGM App can collect, analyze, transform and display the current data collected by the transmitter, store the collected data and generate reports.

Users can get the latest version of User Manual of POCTech CGM App at www.poctechcorp.com

Users can contact the supplier to obtain POCTech Care software for Windows, which is used to input the stored CGM data into the computer through the USB cable to display the glucose data collected by the CGM.

2.5. Receiver Menu Description (Fig. 3):


Figure 3. Receiver Display

📶 Sensor ON/OFF:

To turn transmitter data communication link ON/OFF (Fig.4): In normal operation mode, the “Sensor On” item is not necessary. The sensor/transmitter connection initiates via the “Scan to search” operation. When a sensor is in data link status, the “Sensor On” item is not operable. The “Sensor Off” is available only for particular cases when the monitoring needs to be interrupted, or stopped before the end of the nominal 14-day period.

The “Sensor Signal” is for sensor status check purposes (trouble shooting) when abnormality is suspected. The Sensor Menu can be opened by long pressing the icon (more than 3 seconds).

Figure 4. View Signal Screen

⚙️ Parameters Settings (Fig. 5):

Date/Time;
Keypad Tone;
Alert On/Off (Audio/Vibrator/Audio+Vibrator);
Alert Setting (Hyperglycemia/Hypoglycemia threshold values);
Language.


Figure 5. Parameters Settings Menu

 **Events Input(Fig. 6):**
 RefBG: User enters Blood glucose (BG)
 Event Input: User Logs:
 Exercise;
 Meal;
 Medication;
 Insulin.


Figure 6. Events Input

 **Data Review (Fig.7):**
 Review Data: Review historical sensor data.
 24 Hours Data: Display a 24-hour Overlay of the selected sensor data file.


Figure 7. Data Menu

Section 3. Preparation before Use

3.1. Preparation before Use

The Sensor:

Check the sterile packaging and the expiration date to make sure that the sensor is in good condition.

The Transmitter:

Make sure that one new CR1620 (3V) cell button battery is available for each sensor when beginning the monitoring. Prepare the transmitter (Fig.8): Install a CR1620 battery in the transmitter.


Figure 8. Install Transmitter Battery

Make sure the battery's positive side is upward.
Tighten the battery chamber lid with a coin (Fig.9).


Figure 9. Tighten the Battery Lid

The transmitter starts wireless broadcasting. Its ID number will be displayed on the nearby receiver (When the receiver is in Search Mode).

Figure 10. Transmitter ID Number

The Receiver:

Press and hold the power switch for more than 3 seconds to turn the receiver power ON. Verify that the date and time are correct, and that the battery indicator indicates sufficient battery power.

When the transmitter has been successfully started, its ID number will be displayed on the receiver. Verify the number of the transmitter. Touch the number on the screen to select the transmitter. And then “Yes” to connect the transmitter. The Bluetooth antenna icon appears on the receiver screen, indicating the data communication has been established. Then choose “Sync” to synchronize the transmitter and the receiver (Fig. 11).


Figure 11. Identify Transmitter and Connect

Alcohol Pads and Medical Adhesive Dressing:

Have sufficient supplies available for skin cleaning, disinfection, and transmitter protection.

Section 4. Product Operation and Notice

4.1. Using a New Sensor

Prepare the Receiver:

When the receiver is turned ON, it searches for all nearby transmitters and displays a list of transmitters (or other Bluetooth devices nearby).

Check receiver battery. Charge the battery if necessary.

Set Date and Time, Alerts Modes, Hyper/Hypo Thresholds, and Language

Prepare the Skin:

Choose one upper arm or middle abdomen for sensor insertion (Fig.12). Avoid scars or abnormal areas.

Use alcohol swabs (or pads) to clean and disinfect the skin.


Figure 12. Site Selection for Monitoring

Prepare the Sensor:

Caution

Examine the sensor sterile packaging.
Do not use if the package has been damaged.

Caution

Always check the expiration date on the sensor packaging before opening the sensor. Do not use expired sensors. Tear open the sterile package to take out the sensor assembly. Remove the protection film on the sensor base adhesive tape (Fig.13).


Figure 13. Peeling off Protection Film

Apply the Sensor:

Place the sensor on the cleaned skin.

Press the sensor base and the adhesive tape to achieve a firm adhesion(Fig.14).




Figure14. Press the sensor plate firmly to ensure adhesion

Remove the protection cap on the release button (Fig 15).

Remove the protection cap on the release button (Fig 15)

Bend down the inserter handle until the front sensor base plate pulls the skin up. (Figs. 16a).

Slika 16a.

Push the release button to trigger the insertion mechanism. The sensor should be inserted in the manner shown in Figure 16b.


Figure 16b. The Correct Insertion Angle

After sensor insertion, pull back to remove the inserter lock pin with the thumb of hand holding the inserter (Fig 17).


Figure 17. Remove the Lock Pin

Then, press the arms on the sides of the inserter to disengage the sensor. Remove and discard the inserter in medical sharp disposal (Fig. 18).


Figure 18. Disengage and Separate the Inserter

Press down the Sensor Connection Block (Fig. 19):
Push the sensor connection block down to lock in the base.

Figure 19. Push the Sensor Connection Block into Position

4.2. Initialize the Sensor

The screen then displays “Start Sensor? Yes/No”. When choosing “Yes”, the system enters next Initialization Mode and begins to store sensor data (Fig.20).


Figure 20. Starting Sensor Data Acquisition

Install the transmitter to the sensor: Slide the transmitter onto the sensor base and lock (Fig. 21).


Figure 21. Sliding Transmitter onto Sensor and Lock

Note

Place a medical adhesive dressing tape (3x4”, water proof) to cover the entire transmitter/sensor assembly if necessary.

4.3. Sensor Initialization and First Ref BG

The sensor requires a 1-hour initialization period. The countdown is displayed on the screen.

Caution

Due to physiological differences, blood glucose values within 2 hours after a meal or a high sugar drink cannot be used as calibration references. It may result in significant error!

Note

It is very important that you plan the time to begin your 1-hour initialization. The best time to implant the sensor is immediately after breakfast or lunch so that the first calibration (after a 1-hour initialization) will be before next meal when the blood glucose is in a steady state.

Once the initialization is complete, a reference finger blood glucose (BG) is required. Take a finger stick blood glucose (BG) measurement, Select “”, and then “Ref BG Input” to enter the finger BG (Fig.22).

Figure 22. Entering Reference Blood Glucose

The sensor glucose value and a continuous glucose graph will be displayed after the BG entry (Fig. 23).


Figure 23. Receiver Display for Normal Glucose Monitoring

The sensor requires reference finger blood glucose (BG) for calibration. Follow the next Calibration Table to measure BG and enter to the receiver within 5 minutes.

	1 st	2 nd	Remark
1st Day	First BG	6 – 10 hrs after First BG.	The second BG should be before dinner or before bed.
2nd Day	a.m. Fasting	Before Dinner	
Days 3 – 14	a.m. Fasting	NA	One BG daily for days 3 – 14

Note

You must enter the fasting finger BG before breakfast. Do not enter BG readings within 2 hours after a meal.

4.4. Events Input:

The user should maintain a full activity log by entering all meals, exercise, medicine, and insulin, during the entire sensor wearing period (Fig.24).


Figure 24. Events Log Menu Items

4.5. Features of the Receiver Touch-Screen Display

The Default Display: The default display includes a 5-hour glucose graph and a glucose value. The unit for glucose concentration is mg/dL (mmol/L is available in countries where this unit is used). Normal graph and glucose concentration are displayed in white. Press the HOME key on any screen to return to the default display.

Arrows: When an up or down arrow appears next to the glucose value, it indicates the glucose trend direction of change.

Hypoglycemia: The graph turns to red accompanied with an alert (beep or vibration) when the sensor glucose crosses the pre-set threshold.

Hyperglycemia: The graph turns to yellow accompanied with an alert (beep or vibration) when the sensor glucose crosses the pre-set threshold.

Events Marks: When the user enters an event log (👤), the corresponding event mark appears on the glucose graph (Fig.24).

Zoom-in or Zoom-out the timescale: Touch the bottom of the screen with a finger each time will cause the time scale to zoom out 1/2; Touch the top of the screen with a finger each time will cause the time scale to zoom in doubles.

Left/Right Turn Pages: Touch the left or right sides of the screen moves one page each time.

Green colored Date/Time Means Review Mode: When you review the graph, the date/time display turns green, indicating that the time is not the current clock time.

A Green vertical line and green glucose value is for review point value: Long press on any point of the glucose graph will display a green line. The glucose value of this point is displayed in green.

24-hour Overlay: Touch the data icon (📊) to open the data review menu. Select the “24 hours data” to open a list of sensor data files. Select the last file for the current sensor. A 24-hour overlay plot will be displayed.

Alert On and Off: The beep and vibration alerts are provided to you for unattended alert or reminder when an abnormal situation occurs. To set the Alert Off, use the “Settings” menu (⚙️) to select Alert On/Off, then choose “off”, “OK”, then “confirm OFF”.

4.6. Care and Maintenance during Monitoring

Read the User Manual before using and operating the system. Users including physicians, nurses, and patients must be trained.

Keep the receiver within 4 meters. The designated effective distance of the wireless communication is 4 meters. It is recommended that you carry the receiver at all times. The alert function will not be available when the transmitter is not linked.

Monitor the receiver battery indicator. Charge the battery whenever needed. Use the USB cable provided with this system to charge the receiver. When the instrument has not been used for over a month, it should be charged once. Keep away from fire and high electromagnetic fields. Avoid strong impact. When connecting to a power adapter or a computer, it should be with an EN 62368-1 certified ME system.

Showering: The sensor/transmitter assembly, when properly installed and maintained, is water resistant to an Ingress Protection Rating of IP27 (submerging in water 1 meter/3 feet deep for up to 30 minutes). This means that you can take showers without extra protection. However, due to differences in wearing times and differences in people’s skin conditions, it is strongly recommended that you take necessary precaution and measures when engaging in activities involving water.

Sweating: Excessive sweating should be avoided. It may cause short circuit between sensor connections and adhesive tape to become loose.

Accidental Loss of Sensor: In case the sensor comes off by accident, do not try to re-use. Report the incident or hand it to your doctor.

Outdoor Activities: Avoid vigorous sports activities that might involve physical bumps.

Itchiness: Adhesive dressing tape usually may cause skin itchiness. This can be relieved by changing the cover dressing tape every few days. When changing tape, be careful not to disturb the sensor. If swelling or redness accompanying extreme itchiness happens under the sensor area, the sensor may need to be removed.

Avoiding Extreme Temperature: Try not to enter or stay in extreme temperature that is outside the specified range of the system operation.

Strong Electromagnetic Fields: Avoid strong electromagnetic fields. CT-Scan, X-Ray or MRI: If you must be examined with CT, X-ray, or MRI imaging, remove the transmitter and sensor. Start a new sensor after the examination.

Effect on other Medical Equipment: Under extreme circumstances, nearby medical equipment may be affected. Refer to IEC declaration.

Abnormal Signal: System malfunction, short circuit, sweating, water ingress, etc., may cause signal abnormality. If the situation persists, remove the sensor and report the incident.

The receiver should not be placed close to other similar systems to avoid interference. Insulin pump users should not place the sensor on the same side of the abdomen as the insulin infusion line.

The data system only functions with the designated sensors. Refer to related announcements by the manufacturer regarding sensor compatibility in the future. Do not connect sensors with other unapproved devices or networks.

The information displayed by the receiver is intended to supplement instead of replacing a registered blood glucose measurement device or instrument. It is intended to provide information on continuous changing glucose concentration and changing trends.

Avoid contact by small children. The system contains small parts that may present a choking risk if swallowed.

4.7. Ending Your Continuous Glucose Monitoring

The system is programmed to collect data for up to 14 full days. Data collection automatically terminates when the 14-day acquisition is complete.

Turn off the sensor from the Sensor Menu “” and remove the sensor together with the transmitter. When removing the sensor from the skin, peel and lift the sensor base adhesive from one edge and pull up it firmly until the entire sensor/transmitter assembly is separated from the skin. Do not discard the transmitter.

Check the skin site carefully for bleeding, infection, inflammation etc., if there is a concern, consult medical personnel.

Continue to use your BG meter for self-management of diabetes. Symptoms related to low or high blood glucose levels should not be ignored. If you have symptoms of low or high glucose, use your blood glucose meter to verify.

4.8. Waste Disposal

Sensor: Dispose according to medical waste regulations.

Transmitter: Can be cleaned and reused. May be disposed as medical electronics.

Receiver: Can be cleaned and reused. May be disposed as medical electronics.

Battery: Disposal of the battery should meet local law and regulatory requirements.

Section 5. Troubleshooting and Maintenance

5.1. Troubleshooting

Receiver Memory Over-Flow. The receiver's memory stores no less than 10 cycles of sensor' data. However, it is recommended that you upload the sensor data to your computer (the desktop software can be obtained via customer service or dealers).

Receiver automatically checks the memory space when it starts a new sensor. If there is not enough memory for a sensor's full data it will prompt user to delete old data files to free up memory space.

To delete a data file, go to data icon< (I,II). Select "Review Data" to display a list of data files. Long press the file name you want to delete and respond by selecting "Yes" to confirm the deletion.

Abnormal Sensor Signal. When the sensor signal is out of the expected range, the receiver will sound an alert. The glucose display may be turned off. This usually occurs when water causes a short circuit. A partial pulling off of the sensor may be another reason for alert for abnormal signal.

Sensor Failure. When a sensor failure has been confirmed by the system, the receiver will display "Sensor failure. Please end the session."

Issue establishing initial connection. If there is an issue establishing initial connection between the transmitter and receiver, try powering the receiver off then on. If the issue still exists, replace the transmitter battery.

Transmitter Low Battery. The receiver receives the transmitter battery information when establishing initial connection. A failure message will be displayed and a connection cannot be established if the transmitter battery is lower than 2.9V. Simply replace a CR1620 battery and reconnect.

Receiver Low Battery. The receiver has a rechargeable battery. The battery indicator on the upper right corner changes color when the battery is low. Make sure to check it regularly and charge the battery as needed. The battery normally has a remaining operational time of more than 10 hours even when the indicator color turns orange, allowing enough time to find a chance to charge.

Note

When charging the receiver, a green battery indicator indicates the charging is in progress. It turns white when the charging is complete. DO NOT unplug until the battery indicator turns WHITE.

Hyperglycemia Alert. The high glucose Alert is triggered when the sensor glucose is over the pre-set threshold value. It is advisable not to set the hyper threshold below 200 mg/dl (11.1mmol/L) to avoid frequent alerts. It is strongly recommended that a physician's advice is sought.

Hypoglycemia Alert. Hypoglycemia Alert is the most important feature of the CGM device for patients who may have hypoglycemic unawareness. The low glucose Alert is triggered when the sensor glucose is lower than the hypoglycemia threshold (the default setting for hypoglycemia is 70mg/dL, or 3.9mmol/L). However, the registration condition of this device requires a finger BG measurement to confirm the situation when a hypoglycemia event is detected by the sensor.

Blue-Tooth Communication Interruption. If the receiver alerts for lost communication during normal monitoring, it is generally due to an out-of-range condition between the receiver and the transmitter. Bring them closer together to re-connect.

Unexpected Sensor Removal. Although unlikely, the sensor may accidentally fall off, usually due to excessive sweating or movement, or for people who have oily skin that makes it difficult for the sensor adhesive to adhere firmly. Do not re-use the sensor. Start a new sensor.

Sensor Glucose Obviously Different From Experienced Judgment. When the sensor reading appears erroneous based on personal experience, a finger stick BG measurement is recommended.

Water Ingression. See abnormal sensor signal.

Receiver Time Difference. You may notice that the time on the receiver is different from that of the actual clock. It is normal as long as the time difference is within a few minutes. You are not allowed to reset the time during a sensor monitoring session.

Damage or wear to the Transmitter Battery Lid. An additional battery lid is provided in the transmitter package. It should be used whenever deformation in the O-ring or damage to the lid is observed.

Other Problems. Contact customer service.

Error Code	Receiver Display	Reason	Measure
E01	Receiver battery low! Please charge the battery!	The receiver battery power is too low.	Recharge the receiver.
E02	Memory low, Please delete data files!	Receiver automatically checks the memory space when it starts a new sensor. If there is not enough memory for a sensor's full data it will prompt user to delete old data files to free up memory space.	To delete a data file, go to data icon (📁). Select "Review Data" to display a list of data files. Long press the file name you want to delete and respond by selecting "Yes" to confirm the deletion.
E03	Transmitter battery low. Please replace!	The receiver receives the transmitter battery information when establishing initial connection. A failure message will be displayed and a connection cannot be established if the transmitter battery is lower than 2.9V.	Simply replace a CR1620 battery and reconnect.
E04	Search for transmitter failed. Please Restart the receiver!	Blue Tooth communication error. Battery issue.	If there is an issue establishing initial connection between the transmitter and receiver, try powering the receiver off then on. If the issue still exists, replace the transmitter battery.
E05	Data link lost. Recovering...	Out of range. Strong environmental interference.	Bring them closer together to re-connect. Restart receiver.
E06	Sensor Signal Error! Recovering...	When the sensor signal is out of the expected range, the receiver will sound an alert. The glucose display may be turned off. This usually occurs when water causes a short circuit. A partial pulling off of the sensor may be another reason for alert for abnormal signal	Wait for recovering.(Also reference E08)
E07	Blood glucose is unstable. Please recalibrate when glucose becomes stable.	Blood glucose change in the last 15 minutes is too fast. The measured finger glucose may not be reliable.	Please recalibrate when glucose becomes stable. Do not take a finger blood glucose measurement when the receiver is showing a rapid glucose change in the last 15 minutes. Wait until blood glucose becomes stable to take next measurement for calibration.
E08	Sensor failure. Please end the session.	The anomaly change of current persists over 15 minutes.	Turn off the sensor from the Sensor Menu and remove the sensor together with the transmitter.
E09	Transmitter failure. Please end the session.	Transmitter damage.	Please end the session.

5.2. General Maintenance of the Device - Inspection and Maintenance

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No modification of this equipment is allowed.

Only the manufacturer should examine and provide parts and repair.

Maintenance and repair during use is prohibited.

Transmitter batteries are one-time cell button batteries. They may be changed by the user. The battery should not be installed when not in use. The receiver battery is a rechargeable lithium battery. Changing of the receiver battery, if necessary, must be performed by designated representatives.

Section 6. Specifications and EMC Declaration

6.1. Specifications

- Test Sample: The test sample is the subcutaneous tissue.
- Product Parameters:
 - i. Sensor Shelf life: 12 months;
 - ii. Glucose Range: 1.7-25mmol/L(30 - 450mg/dL);
 - iii. Lab Test Accuracy:

iv. Range	Specification
< 4,2 mmol/L (< 75 mg/dL)	MAD*≤0.83mmol/L (≤15 mg/dL)
≥ 4,2 mmol/L (≥ 75 mg/dL)	MARD**≤12%

*MAD, Mean Absolute difference;

**MARD, Mean Absolute Relative Difference

Clinical data

A clinical trial for continuous glucose monitoring system (CGMS) was completed in 2015 by POCTech, involving 80 patients in two hospitals. An overall MARD (vs finger sticks) of 9.61% was achieved. The sensor also showed an overall MARD of 8.67% (vs venous blood glucose).

- v. Transmitter-Receiver Communication Range: 4 meters (~12 feet);
- vi. Frequency: 2.4 GHz; Transmitter E.I.R.P: 1.3 dBm: Receiver E.I.R.P: 1.8dBm
- vii. Transmitter Weight: 7.6 grams (~0.27ounces);
- viii. Receiver Weight: 98 grams (~3.46 ounces);

- Product Environment Requirement:

i. Transportation Conditions

Condition	Data System	Sensor
Temperature/°C Humidity /RH, %	-15 – 45 °C ≤ 93 % RH	2 – 45 °C, ≤30 days

ii. Storage Conditions

Condition	Data System	Sensor
Temperature/°C	-15 – 45 °C	2 – 30 °C
Humidity /RH, %	≤ 93 % RH	NA
Electrostatic	Protected	NA

iii. Operating Condition

Condition	Data System	Sensor
Transmitter Battery	DC 3,0 V, (CR1620)	NA
Receiver Battery	DC 3,7 V	NA
Pressure/kPa	70 – 106 kPa	70 – 106kPa
Temperature/ °C	5 – 40 °C	15 – 40 °C
Humidity(RU)/%	≤ 93 % RH	NA

6.2. EMC Declaration

Guidance and Manufacturer's Declaration – Electromagnetic emission		
CT14 Continuous Glucose Monitoring System is expected to be used in the following environment. Buyers and users should ensure that such electronic conditions are met.		
Emission Test	Conformity	Electromagnetic Environment – Guidance
RF emission CISPR11	Group 1	CT14 system adopts emission only for its internal communication. The energy level is very low and is not expected to cause any interference to surrounding electronic equipment.
RF emission CISPR11	Type B	/

Guidance and Manufacturer's Declaration -Electromagnetic Immunity

CT14 Continuous Glucose Monitoring System is expected to be used in the following environment. Buyers and users should ensure that such electronic conditions are met.

Immunity Test	IEC60601Electrical Level Test	Met Electrical Level	Electromagnetic Environment – Guidance
Electrostatic discharge IEC61000-4-2	±8KV contact discharge ±15KV air discharge	±8KVcontact discharge ±15KVair discharge	The ground surface should be wood, concrete, or ceramic. If the ground is of synthetic materials, the relative humidity should be at least 30%.
Power frequency magnetic field (50Hz,60H) IEC61000-4-8	3 A/m	30 A/m	Power frequency magnetic field should be of the typical levels of business or hospital environment.
RF emission IEC61000-4-3	3 V/m, 10 V/m 80MHz-2,7GHz	3 V/m, 10 V/m 80MHz-2,7GHz	Portable and mobile RF communication equipment (including any part of the CT14 system) should be no closer than that of recommendation below. Recommended distance calculation: $d = 1,167 \times \sqrt{P}$ 80 MHz – 800 MHz $d = 2,333 \times \sqrt{P}$ 800 MHz – 2,7 GHz P—The maximum output power level of the Transmitter based on the manufacturer, unit (W). d—Recommended distance, Unit (m). The field intensity of an immobilized RF Transmitter is determined by the a factor. It should be below the electrical level at each b level. Interference may occur near equipment with such label. (())

Note 1: At 80MHz and 800MHz frequencies, use higher frequency formula.

Note 2: The guidance may not be suitable for all circumstances. Electromagnetic transmission is affected by buildings, objects and human bodies.

- Fixed position Transmitters, such as cordless phones, hubs for wireless networks, amateur shortwave radios, FM,/AM radios, and TVs, etc., emit electronic fields whose intensity may not be accurately known. To evaluate the intensity of such fields, observation should be done regarding such environment. If higher than the above tested electric levels are observed, and equipment interferences confirmed, supplemental measures may be necessary, such as adjusting direction or location.
- The field intensity should be below [V1]V/m in the entire range.

Recommended Distances for Portable or Mobile Communication Equipment and Systems-Non-Life Supporting Equipment and Systems.

Recommended Distances for Portable RF Equipment with CT14 system

CT14 Continuous Glucose Monitoring System is expected to be used in an environment where RF interferences are controlled. Based on the maximum output power levels, buyers or users may use the below calculation to estimate the minimum distances from an RF emitting device required to use the CT14 system.

The maximum nominal output power level (W)	$d = 1,167 \times \sqrt{P}$ (80 MHz – 800 MHz)	$d = 2,333 \times \sqrt{P}$ (800 MHz – 2,7 GHz)
0,01	0,12	0,23
0,1	0,37	0,74
1	1,17	2,33
10	3,69	7,38
100	11,67	23,33

For equipment whose output levels not listed in the above table, an estimate can be made using the equation, where the unit for d is meter (m), and p (the maximum power) is watt (W).

Note 1: At 80MHz and 800MHz frequencies, use higher frequency formula.

Note 2: The guidance may not be suitable for all circumstances. Electromagnetic transmission is affected by buildings, objects and human bodies.

Basic properties tested in the electromagnetic tests:
DC current measurement in the range of 0~50nA, within±1nA.
DC current measurement in the range of 50~370nA, within ≤2.0%.

6.3. EU Radio Equipment Directive (RED)Declaration

Hereby, POCTech Co.,LTD declares that the radio equipment type CT-100C14 Transmitter and CT-100BD Receiver is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address:www.poctechcorp.com/CGMS/declarationofconformity

Section 7. Warranty

7.1. Contacts and Numbers:

- Customer Service

Customer Service	Zhejiang POCTech Co., Ltd.
Phone	86-400-118-8528
Web	www.poctechcorp.com

7.2. Warranty

Warranty period for the transmitter and receiver: one year from the date of purchase. The following conditions may incur charges even within one year:

- Damage or malfunction caused by false use or misuse;
- Damage or malfunction caused by force including but not limited to: fire and earthquake.
- Damage or malfunction caused by moving or shipping after purchase.
- Damage or malfunction caused by other causes not attributed to the manufacturing.
- POCTech Co., Ltd. reserves right of explanation to the above terms and conditions.

Zhejiang POCTech Co., Ltd. No. 1633
Hongfeng Road, Building 11 & 12
Huzhou City 313000 , Zhejiang ,China

Prolinx GmbH
Brehmstr. 56, 40239 Duesseldorf ,
Germany

Version:V20211230

Zhejiang POCTech Co., Ltd.
No. 1633 Hongfeng Road, Building 11 & 12,
Huzhou City 313000 ,
Zhejiang ,China
86-400-118-8528